

TRYP- Complementation

	694/A74	694/A10	695/A28	1150/A62	1148/A38	1147/A20	1149/A50	1152/K142	1151/K976
1) 694/A74	1	2	3	4	5	6	7	8	9
2) 694/A10	responses to antibiotic test		11	12	13	14	15	16	17
3) 695/A28			18	19	20	21	22	23	24
4) 1150/A62				25	26	27	28	29	30
5) 1148/A38					31	32	33	34	35
6) 1147/A20						36	37	38	39
7) 1149/A50							40	41	42
8) 1152/K142								43	44
(1151) 1151/K976									45

responses to
antibiotic test

(OK -)
(OK comment OK)

replaced

1149?
1151?

Test was run on following medium - with/without NZ Case. Results were best on medium containing NZ Case. Therefore, use following for Complementation:

- Distilled water 1000 ml
- Bojels (50X) N medium 20 ml
- glycerol 1.0 ml
- glucose 4.0 gm
- yeast 20.0 gm
- Sheffield's NZ Case 0.2 gm
- Sorbitol 15.0 gm