

MEDIA FOR SORDARIA FIMICOLA

	CORNMEAL AGAR DEHYDR.	SUCROSE	GLUCOSE	KH ₂ PO ₄	Na-Ac	YEAST EXT.	H ₂ O (CC)	AGAR
CM ⁺ (GENOTYPE DETERMINATION, STOCK CULTURE)	17g	-	2g	-	-	1g	1,000	-
CULTURE (ALL CROSSES)	17g	10g	7g	0.1g	-	1g	1,000	-
DISSECTING (NEW) (GOOD FOR ANY DISSECTION, ESSENTIAL FOR <u>h</u> SPORES)	2g	15g	10g	-	1g	-	500	15g
DISSECTING (TRADITIONAL) (GOOD FOR CONTAMINANT REMOVAL)	8.5g	-	1g	-	3.5g	0.5g	500	15g
TESTING	17g	15g	12g	0.1g	-	2g	1,000	-
HYPERTONIC MOUNTING SOLUTION	-	10g	5g	-	0.7g	-	100	-