

keep stocks at 4 degrees C and passage them every 3 months. We grow liquid and plate cultures for 10-14 days at room temperature (growth is usually obvious by 1 week and plates grown longer than 2 weeks generally dry out). We usually inoculate small agar plates or 30ml liquid vials with 1mL of liquid culture or one colony from a previous plate. Plates can be checked easily for healthy growth because motile zoospores and characteristic sporangia are visible even under low magnification. We also have a protocol for cryopreserving if you would like it. Let me know if you need any additional information.

(use back of page if necessary)

YOUR NAME Erica Rosenblum

DATE 4/16/07

11501  
Erica Bree Rosenblum

email: rosenblum@berkeley.edu

phone: 510-508-2622

Marlby address:

One Cyclotron Road

Bldg 84, room 355

Berkeley CA 94720

up hand  
2007-04-16