

Hi Kevin,

Here are the strains we are sending. The # is meaningless for you, just our cataloging system. Ve+ = veA+ (most of the strains you have at FGSC are veA1).

*Remarks as from Larry Yager*

Let me know if you have questions.

#	Strain	Aux	Status	spore	growth	genotype	Remarks
1	WIM001	mm	+	green	normal	aco-3; ve+	OK; temperature-sensitive (ts)
25	WIM 113	Paba	+	green	Sick	paba; acoA49; ve+	(ts) spores look yellow
2	WIM 002	mm	+	green	normal	aco-14; ve+	spores look yellow (ts)
26	WIM 126	Paba	+	yellow	normal	yA2; pabaA1; ve+	KEEP THIS ONE - IT RESPONDS TO BLUE LIGHT while FGSC 4 DOES NOT.
3	WIM 003	mm	+	green	sick	acoA49; ve+	Also a KEEPER. The original competence mutation. Looks yellow; it's green. Grows really bad, tiny colonies. Try growing at 28 C. (ts)
4	WIM 004	mm	+	green	normal	aco-55; ve+	OK (ts)