

A2066 MG151: *pyrG^{AF}* insertion (*pyrG89*); *nimT23*; Δ *Ku70::argB* (*argB2*); *pyroA4*; *pabaA1*; *sE15*; *wA3*.

A2067 MG411 Δ *An-swd1::pyro^{AF}* (*pyroA4*); *H3K4R::pyrGAF* (*pyrG89*); *sE15*; *nirA14*; *wA3*; *fwA1*; *chA1*.

Question marks indicate mutations or alleles that may be present in the strain and have not been tested for.

Strains in Red need to be grown on YAG + 1M Sucrose (otherwise they donot conidiate well)