

3961 and 3962 - "(i) Both give am i recombinants on cross to am. Growth on N-free Vogel's liquid minimal (sucrose) with 0.01 M L-proline - 15 ml liquid in 100 ml erlenmeyer. 2½ days - no shaking. 0.05 ml drop of inoculum ca 10⁶ conidia/ml. The wild gives about 15 mg and i about 2 mg. Proline agar plates do not differentiate as well, if at all. (ii) Minimal (Vogel's) agar sorbose plates with 0.2 mM DL-p-fluorophenylalanine. Wild type hardly grows at all. i is a bit inhibited but shows a good colony from a loop on spot inoculum in 2 days. This is very clear. As a final test cross to am. Around 5% of the ascospores are am i doubles and grow extremely badly on all media - you have to have the NH₄⁺ out -- e.g. glutamate as sole N-source -- to get reasonable growth." (J.R.S. Fincham).

3961 and 3962 - "(i) Both give am i recombinants on cross to am. Growth on N-free Vogel's liquid minimal (sucrose) with 0.01 M L-proline - 15 ml liquid in 100 ml erlenmeyer. 2½ days - no shaking. 0.05 ml drop of inoculum ca 10⁶ conidia/ml. The wild gives about 15 mg and i about 2 mg. Proline agar plates do not differentiate as well, if at all. (ii) Minimal (Vogel's) agar sorbose plates with 0.2 mM DL-p-fluorophenylalanine. Wild type hardly grows at all. i is a bit inhibited but shows a good colony from a loop on spot inoculum in 2 days. This is very clear. As a final test cross to am. Around 5% of the ascospores are am i doubles and grow extremely badly on all media - you have to have the NH₄⁺ out -- e.g. glutamate as sole N-source -- to get reasonable growth." (J.R.S. Fincham).