

FGSC #	Lyophilized	Viability	Silica Gel	Viability
4478	2/16/84	4/2/84 OK	2/16/84	4/2/84 OK
	Genotype responds to cit+ urid+lys+nic+ind not singly		Comments	
FGSC #	Lyophilized	Viability	Silica Gel	Viability
4479	2/16/84	4/2/84 OK	2/16/84	4/2/84 OK
	Genotype slowly → responds to: cit+urid+ lys+nic		Comments <u>check (L)</u> un OK - seems to be lys ⁺ and ind ⁺	
FGSC #	Lyophilized	Viability	Silica Gel	Viability
4480	3/15/84	4/2/84 OK	3/15/84	4/2/84 OK
	Genotype cit - urid - ind - all 3 +		Comments	
FGSC #	Lyophilized	Viability	Silica Gel	Viability
4481	3/15/84	3/24/84 OK	3/15/84	3/22/84 OK
	Genotype cit - urid - cit+urid+		Comments	
FGSC #	Lyophilized	Viability	Silica Gel	Viability
4482	2/16/84	3/20/84 OK	2/16/84	3/20/84 OK
	Genotype cit - nic - cit+nic+		Comments	
FGSC #	Lyophilized	Viability	Silica Gel	Viability
4483	3/15/84	3/29/84 OK	3/15/84	3/22/84 OK
	Genotype cit - nic - urid - all 3 +		Comments	
FGSC #	Lyophilized	Viability	Silica Gel	Viability
4484	3/15/84	3/22/84 OK	3/15/84	3/22/84 OK
	Genotype min - ind +		Comments	